



MIRION
TECHNOLOGIES

DETECTORS

NaI Detector Accessories

7419 Shield/Collimator



FEATURES

- All attitude
- Compact: swivel legs and folding handle
- Removable collimator: 44 mm (1.75 in.)
- Removable tin filter: 1 mm thick (0.040 in.)
- Lead shielding: 8 mm wall x 122 mm long (0.32 x 4.8 in.)
7419BS: 8 mm wall x 76 mm long (0.32 x 3.0 in.)
- Shipping weight: 6.5 kg (14 lb)
- Five models – specify when ordering:
 - 7419A for 802-2 x .5
 - 7419B for NAID
 - 7419BS for NAID (with detachable legs)
 - 7419C for 802-2 x 2
 - 7419E for IPROS-2

ISO 9001 MERIDEN CT
SYSTEM SMYRNA GA
CERTIFIED OAK RIDGE TN
CONCORD ON

M1184 – 08/2012

Copyright © 2019 Mirion Technologies, Inc. or its affiliates. All rights reserved. Mirion, the Mirion logo, and other trade names of Mirion products listed herein are registered trademarks or trademarks of Mirion Technologies, Inc. or its affiliates in the United States and other countries. Third party trademarks mentioned are the property of their respective owners.



MIRION
TECHNOLOGIES

www.mirion.com